



MEASURE THE ANS  
CONTROL NOCICEPTION  
IMPROVE OUTCOMES

## ANI in the ICU

Analgesia  
Nociception Index



## The problem

In the intensive care unit (ICU) there are patients requiring mechanical ventilation and sedation. Analgesia is key to help patients with a high respiratory drive adapt to the ventilator and prevent patient-ventilator asynchrony (PVA).

To achieve good patient-ventilator synchrony optimal balance between sedation and analgesia is required. This balance can be achieved more effectively by monitoring the effect of each procedure with:

### An EEG monitor

To control the desired level of consciousness.

### The ANI monitor

To detect noxious stimulation. ANI may help to guide analgesia [8].



ANI MOC-9 Module



ANI-MR Module



ANI Monitor V2



Several publications have proven that, in mechanically ventilated patients, an analgosedation protocol will improve pain scores, **reduce sedatives (54%)** and their overall consequences as delirium, **decrease length of stay (50.8 h)**, **decrease time on ventilation (45.5h)** and decrease associated medication costs [1]. However, overuse of analgesia, especially opioids, must also be avoided as 20% to 70% of patients that receive opioids will have side effects such as; respiratory depression (27%), gastrointestinal problems (41%), delayed emergence (28%) and tachyphylaxis (37%) [2]. Hence, it is key to be able to specifically titrate the analgesic requirements of each patient.

### The benefits of an analgosedation protocol

**54%**

reduction  
of sedatives


**50,8h**

decrease length  
of stay

**45,5h**

decrease time  
on ventilation

### % of opioid-overdosing side effects




# Our solution

## Measuring the autonomic nervous system (ANS) with the analgesia nociception index (ANI)

Heart rate variability (HRV) is a well-known phenomenon that controls the regulation of the cardiovascular system via the ANS. By analyzing the HRV oscillations it is possible to measure the activity of the sympathetic and parasympathetic nervous systems [3].

All mammals exhibit what is known as respiratory sinus arrhythmia. This process is controlled by the sympathetic and the parasympathetic branch of the nucleus accumbens and nucleus ambiguus, causing the heart to beat slower during expiration, due to vagal activation, and faster in inspiration due to sympathetic innervation [4].

### How to obtain ANI from ECG




\*AUC : calculation of area under the curve of the respiratory pattern

### The ANI value measures the relative parasympathetic activity of the patient

The total energy of the ANS is shown in the screen of the monitor. During general anesthesia, an ANI range between 50-70 relates to adequate analgesia, meaning that antinociception is adequate and that parasympathetic activity is mildly predominant over sympathetic activity. When the ANI value falls below 50, the occurrence of a hemodynamic response within the following 5 minutes is very likely [5]. This information can be used to predict and avoid a hemodynamic response by increasing analgesia.

## Interpretation of ANI




nb : the energy must be between 0,05 and 2,5 for an interpretable ANI.

# The result

**ANI can detect nociception [6] and may help to take the decisions to control the autonomic response improving patient outcomes [17]**

In Anesthesiology 2017, Funcke et al. show that ANI is the best in its class for detecting nociception versus other parameters, especially under opioid conditions [6]. While titrating analgesics based on heart rate and blood pressure might be common, Funcke et al. demonstrated that ANI is more sensible and specific than haemodynamics to detect nociception [6].

Validation of Innovative Techniques for Monitoring Nociception during General Anesthesia: A Clinical Study




Picture modified from Funcke et al. [6]




Several publications indicate ANI can detect nociception and be safely used in the ICU setting for guiding analgesia in deeply sedated patients, even under catecholamine treatment [10-12].

ANI and BIS variations in supine and prone position during closed-tracheal suction in sedated and myorelaxed ICU patients with severe COVID-19: A retrospective study.



Picture modified from Boselli E, et al. [11]

Measurement of Heart Rate Variability to assess pain in sedated critically ill patients: A prospective observational study.



Picture modified from Broucqsaault-Dédrie C, De Jonckheere J, Jeanne M, Nseir S. [10]

ANI can detect nociception and correlates with Numerical Rating Scores in postsurgical patients who underwent anesthesia with halogenated agents [13].

## What medical societies say about ANI in ICU

ANI has been pointed out as an interesting technology for nociception detection in the 2018 Society of Critical Care Medicine clinical guidelines for sedation in ICU [14], in the sedoanalgesia and ventilation protocol for Covid-19 patients from Hospital Puerta de Hierro de Madrid [15] and by the SEMICYUC (Spanish Society of Intensive Care Medicine) [16]. MDoloris technology is being increasingly used worldwide for the detection and management of nociception.

To find out more about MDoloris technologies please visit our website <https://mdoloris.com/>



## The main benefits of using ANI technology

ANI shows the sympathetic/parasympathetic allowing taking better decisions to control nociception improving outcomes.



Predictivity of hemodynamic reactivity [5]



Helpful to diagnose the etiology of a haemodynamic event [7]



Predict post-extubation pain [9]



Refine opioids titration [8]

## Testimonials



*“ ANI monitor is very useful for monitoring the patient’s autonomic nervous system, especially in those undergoing mechanical ventilation, to titrate analgosedation. This allows good patient synchrony to the ventilator, avoiding the side effects of overdose, saving drugs, and allowing for earlier weaning from mechanical ventilation. ”*

**Dr. Carlos Chamorro Jambrina**

**Intensive Care Unit -  
Hospital Puerta de Hierro Majadahonda, Madrid**



*“ The measurement of parasympathetic tone with ANI to continuously evaluate the level of analgesic comfort of the patient, optimizes the consumption of opiates, achieving adequate analgesia, and with a high negative predictive value that helps to differentiate the etiology of hemodynamic adverse events in the ICU mainly in patients under mechanical ventilation. It allows managing individualized doses to each patient managing to avoid sub or supratherapeutic doses, which has been of particular relevance during the COVID-19 pandemic due to the decrease in the supply of medicines for sedoanalgesia. ”*

**Dr. Angel Augusto Pérez Calatayud**

**Chef of Critical Area Division,  
Hospital General de México**

# Bibliography:

- [1] Andrew C Faust , Pearl Rajan, Lyndsay A Sheperd, Carlos A Alvarez, Phyllis McCorstin, Rebecca L Doebele Impact of an Analgesia-Based Sedation Protocol on Mechanically Ventilated Patients in a Medical Intensive Care Unit *Anesth Analg*, 903-9 Oct 2016
- [2] B V Christensen , L P Thunedborg Use of Sedatives, Analgesics and Neuromuscular Blocking Agents in Danish ICUs 1996/97. A National Survey *Intensive Care Med* 25 (2), 186-91 Feb 1999
- [3] European Heart Journal (1996) 17, 354-381 Guidelines Heart rate variability Standards of measurement, physiological interpretation, and clinical use Task Force of The European Society of Cardiology and The North American Society of Pacing and Electrophysiology (Membership of the Task Force listed in the Appendix)
- [4] Jeanne M, Logier R, De Jonckheere J, Tavernier B. Validation of a graphic measurement of heart rate variability to assess analgesia/nociception balance during general anesthesia, *IEEE Proceedings* 2009.
- [5] Jeanne M, Delecroix M, De Jonckheere J, Keribedj A, Logier R, Tavernier B. Variations of the Analgesia Nociception Index during propofol anesthesia for total knee Replacement. *Clin J Pain* 2014 Dec; 30(12):1084-8.
- [6] Funcke et al. Validation of Innovative Techniques for Monitoring Nociception during General Anesthesia. *Anesthesiology* 2017 Aug;127(2):272-283
- [7] Logier R, et al. Heart rate variability analysis for arterial hypertension etiological diagnosis during surgical procedures under tourniquet. *Conf Proc IEEE Eng Med Biol Soc* 2011;2011:3776-9.
- [8] Daccache G, et al. A targeted remifentanil administration protocol based on the analgesia nociception index during vascular surgery. *Anaesth Crit Care Pain Med* 2017 Aug;36(4):229-232.
- [9] Boselli E, et al. Prediction of immediate postoperative pain using the analgesia/nociception index: a prospective observational study. *Br J Anaesth* 2014 Apr;112(4):715-21.
- [10] Broucqsaault-Dédrie C, De Jonckheere J, Jeanne M, Nseir S. Measurement of Heart Rate Variability to assess pain in sedated critically ill patients: A prospective observational study. *Plos One* 2016; 11(1): e0147720
- [11] Boselli E, Fatah A, Ledochowski S, Allaouchiche B. ANI and BIS variations in supine and prone position during closed-tracheal suction in sedated and myorelaxed ICU patients with severe COVID-19: A retrospective study. *J Clin Monit Comput*. 2020 Nov 6 : 1-7.
- [12] Ali Jendoubi, Ahmed Abbes, Salma Ghedira, and Mohamed Houissa. Pain Measurement in Mechanically Ventilated Patients with Traumatic Brain Injury: Behavioral Pain Tools Versus Analgesia Nociception Index *Indian J Crit Care Med*. 2017 Sep; 21(9): 585-588.
- [13] Abdullayev R, Uludag O, Celik B. Analgesia Nociception Index: assessment of acute postoperative pain. *Rev Bras Anesthesiol*. 2019 Jul - Aug;69(4):396-402. doi: 10.1016/j.bjan.2019.01.003. Epub 2019 Aug 6.
- [14] Devlin JW et al., Executive Summary: Clinical Practice Guidelines for the Prevention and Management of Pain, Agitation/Sedation, Delirium, Immobility, and Sleep Disruption in Adult Patients in the ICU. *Crit Care Med*. 2018 Sep;46(9):1532-1548.
- [15] Carlos Chamorro, Miguel Angel Romera, Sara Alcantara. ICU Department Hospital Universitario Puerta de Hierro. Protocolo de sedoanalgesia para los pacientes con COVID19.
- [16] <https://semicyuc.org/wp-content/uploads/2020/03/PROTOCOLO-ANALGOSED-COVID-SEMICYUC.pdf>
- [17] Henry D. Upton, MBBS, BMedSc (Hons), Guy L. Ludbrook, MBBS, FANZCA, PhD, Andrew Wing, BMBS (Hons), BSci (Hons), FANZCA, and Jamie W. Sleigh, MD Intraoperative "Analgesia Nociception Index"- Guided Fentanyl Administration During Sevoflurane Anesthesia in Lumbar Discectomy and Laminectomy: A Randomized Clinical Trial *Anesthesia & Analgesia* July 2017 doi: 10.1213/ANE.0000000000001984


MDM5  
Mikros Medical Systems

MDM5

11:36:39

Stop

ANI



$\rho\Sigma$

68 m


90 i

Screen shot

ANI navigation

Events

Parameters



Energy 0.35

Quality  
Signal

Good

Reset ECG



ANI



MDoloris Medical Systems SAS  
Biocentre Fleming - Bâtiment C Epi de Soil  
270, rue Salvador Allende  
59120 Loos - France  
Tel: +33 (0)3 62 09 20 80 / Fax +33 (0)9 72 38 75 27  
[www.mdoloris.com](http://www.mdoloris.com)



ANI Monitor V2, ANI MOC-9 and the ANI-MR are class IIa medical device, manufactured by MDoloris Medical Systems. CE evaluation was performed by Bureau Veritas Italy (1370) for the ANI Monitor V2 and the ANI MOC-9. CE evaluation was performed by BSI (2797) for the ANI-MR. © 2021 MDoloris Medical Systems. All rights reserved.  
MD/QUA/EN14.2 v.04